



Ovlašćeni zastupnik i uvoznik:

MEDICALTECH SOLUTIONS DOO
Ćirpanova 15, Novi Sad

Tel. +381 21 2715 767
Mobilni: +381 62 281 640
office@medicaltech.rs
www.medicaltech.rs

Inspiring for Splendid Healthcare

ISO 9001:2015 | ISO 13485:2016 | CE



OSSEOUS PRODUCT CATALOG

- OTOLOGICAL MIDDLE EAR PROSTHESIS
- OTOLOGICAL VENTILATION TUBES
- RHINOLOGY PRODUCTS
- MEDICAL PVA DRESSINGS
- LARYNGOLOGY PRODUCTS
- ENT INSTRUMENTS
- ENT PLASMA SYSTEM - Coaglator II
- ABSORBABLE GELATIN SPONGES

Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2114	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2114B	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2114T	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2114W	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2114BW	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2114BT	1.14	2.4	2.4	2.2	10 Units, ETO Sterile
E2138	1.00	2.4	2.4	2.2	10 Units, ETO Sterile
E2138W	1.00	2.4	2.4	2.2	10 Units, ETO Sterile
E2138T	1.00	2.4	2.4	2.2	10 Units, ETO Sterile
E2138B	1.00	2.4	2.4	2.2	10 Units, ETO Sterile
E2138BW	1.00	2.4	2.4	2.2	10 Units, ETO Sterile
E2138BT	1.00	2.4	2.4	2.2	10 Units, ETO Sterile



Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2137	1.14	2.3	2.3	1.6	10 Units, ETO Sterile
E2137B	1.14	2.3	2.3	1.6	10 Units, ETO Sterile




Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2136	1.27	2.8	2.8	2.2	10 Units, ETO Sterile
E2136B	1.27	2.8	2.8	2.2	10 Units, ETO Sterile



Material : Fluoroplastic (PTFE)

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E2123	1.27	3.0	3.0	2.15	10 Units, ETO Sterile
E2123B	1.27	3.0	3.0	2.15	10 Units, ETO Sterile
E2123T	1.27	3.0	3.0	2.15	10 Units, ETO Sterile
E2125	1.14	2.50	2.50	2.0	10 Units, ETO Sterile
E2125B	1.14	2.50	2.50	2.0	10 Units, ETO Sterile
E2125T	1.14	2.50	2.50	2.0	10 Units, ETO Sterile



Collar Button (ID 1.27mm)
Collar Button (ID 1.14mm)



Donaldson Silicone

Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E3111	1.10	2.3	2.3	2.2	10 Units, ETO Sterile

Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E3126	1.10	3.5	2.5	1.65	10 Units, ETO Sterile



Paparella

Material : Medical Grade Silicone


REF Number	Inner Diameter	Inner Flange Diameter	Total Length	Units/Box
E3115S	1.14	8.0	6.0	10 Units, ETO Sterile
E3115M	1.14	8.0	9.0	10 Units, ETO Sterile
E3115L	1.14	8.0	12.0	10 Units, ETO Sterile



T Tube

Material : Silver Plated Gold

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box
E1170	1.25	2.55	2.55	1.20	10 Units, Gamma Sterile



Gold Vent Tube Type
Collar Button with Wire

B - Blue, T - Tail, W - Wire, All dimensions are in MM

MEDICALTECH SOLUTIONS DOO, +381 62 281 640, office@medicaltech.rs



REF Number	Descriptions	Units/Box
E6149S	External Aluminium Splints Small	10 or 25 Units, ETO Sterile / Non Sterile
E6149M	External Aluminium Splints Medium	10 or 25 Units, ETO Sterile / Non Sterile
E6149L	External Aluminium Splints Large	10 or 25 Units, ETO Sterile / Non Sterile



Material: Medical Grade Silicone

REF Number	Product Name	Specification	Units/Box
E31030	Epistaxis Balloon Catheter Intra-Pack	30 CC / ml	10 Units, ETO Sterile
E31010	Epistaxis Balloon Catheter Post-Pack	10 CC / ml	10 Units, ETO Sterile
E31031	Epistaxis Balloon Catheter Double Balloon	30/10 CC / ml	10 Units, ETO Sterile