

Ovlašćeni zastupnik i uvoznik:

MEDICALTECH SOLUTIONS DOO  
Ćirpanova 15, Novi Sad  
Tel. +381 21 2715 767  
Mobilni: +381 62 281 640  
[office@medicaltech.rs](mailto:office@medicaltech.rs)  
[www.medicaltech.rs](http://www.medicaltech.rs)




# ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY



Graves Vaginal Speculum  
With Fixed or Disposable Smoke Tube

| Fixed Tube Cat # | Disposable Tube Cat # | Size    | Length | Width |
|------------------|-----------------------|---------|--------|-------|
| 130-104          | 130-104 DT            | X-Large | 150 mm | 40 mm |
| 130-100          | 130-100 DT            | Large   | 115 mm | 35 mm |
| 130-102          | 130-102 DT            | Medium  | 95 mm  | 35 mm |
| 130-103          | 130-103 DT            | Small   | 75 mm  | 20 mm |



Disposable Tube



VU-More (Graves) Speculum  
With Fixed or Disposable Smoke Tube

| Fixed Tube Cat # | Disposable Tube Cat # | Size   | Length | Width |
|------------------|-----------------------|--------|--------|-------|
| 130-117          | 130-117 DT            | Large  | 115 mm | 35 mm |
| 130-118          | 130-118 DT            | Medium | 100 mm | 35 mm |
| 130-119          | 130-119 DT            | Small  | 80 mm  | 30 mm |



VU-Max Graves Speculum  
With Fixed or Disposable Smoke Tube

| Fixed Tube Cat # | Disposable Tube Cat # | Size    | Length | Width |
|------------------|-----------------------|---------|--------|-------|
| 130-107          | 130-107 DT            | X-Large | 135 mm | 37 mm |
| 130-106          | 130-106 DT            | Large   | 110 mm | 35 mm |
| 130-105          | 130-105 DT            | Medium  | 100 mm | 35 mm |



Pederson Speculum  
With Fixed or Disposable Smoke Tube

| Fixed Tube Cat # | Disposable Tube Cat # | Size    | Length | Width |
|------------------|-----------------------|---------|--------|-------|
| 130-204          | 130-204 DT            | X-Large | 165 mm | 25 mm |
| 130-200          | 130-200 DT            | Large   | 115 mm | 25 mm |
| 130-202          | 130-202 DT            | Medium  | 100 mm | 22 mm |
| 130-203          | 130-203 DT            | Small   | 75 mm  | 14 mm |



VU-Max Pederson Speculum  
With Fixed or Disposable Smoke Tube

| Fixed Tube Cat # | Disposable Tube Cat # | Size   | Length | Width |
|------------------|-----------------------|--------|--------|-------|
| 130-114          | 130-114 DT            | Large  | 115 mm | 25 mm |
| 130-115          | 130-115 DT            | Medium | 100 mm | 22 mm |



Chinese Pattern Speculum  
With Smoke Tube

| Cat #   | Size   | Length | Width |
|---------|--------|--------|-------|
| 130-255 | Medium | 95 mm  | 28 mm |
| 130-257 | Small  | 70 mm  | 35 mm |

# ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY



Cusco Speculum  
With Smoke Tube

| Cat #   | Size   | Length | Width |
|---------|--------|--------|-------|
| 130-108 | Large  | 95 mm  | 36 mm |
| 130-109 | Medium | 85 mm  | 35 mm |
| 130-110 | Small  | 75 mm  | 32 mm |




Collin Speculum  
With Smoke Tube

| Cat #   | Size   | Length | Width |
|---------|--------|--------|-------|
| 130-111 | Large  | 110 mm | 40 mm |
| 130-112 | Medium | 100 mm | 35 mm |
| 130-113 | Small  | 95 mm  | 30 mm |



Weisman Speculum (Right Open)  
With Smoke Tube

| Cat #   | Size   | Length | Width |
|---------|--------|--------|-------|
| 130-221 | Large  | 115 mm | 35 mm |
| 130-225 | Medium | 95 mm  | 35 mm |



Weisman Speculum (Left Open)  
With Smoke Tube

| Cat #   | Size   | Length | Width |
|---------|--------|--------|-------|
| 130-220 | Large  | 115 mm | 35 mm |
| 130-224 | Medium | 95 mm  | 35 mm |



Weisman Pederson Speculum 90° (Right Open)  
With Smoke Tube

| Cat #      | Size   | Length | Width |
|------------|--------|--------|-------|
| 130-206-90 | Large  | 115 mm | 25 mm |
| 130-207-90 | Medium | 100 mm | 25 mm |



Weisman Pederson Speculum 90° (Left Open)  
With Smoke Tube

| Cat #      | Size   | Length | Width |
|------------|--------|--------|-------|
| 130-208-90 | Large  | 115 mm | 25 mm |
| 130-209-90 | Medium | 100 mm | 25 mm |

# ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY



Weisman Speculum 90° (Right Open)  
With Smoke Tube

| Cat #      | Size    | Length | Width |
|------------|---------|--------|-------|
| 130-227-90 | X-Large | 135 mm | 40 mm |
| 130-221-90 | Large   | 115 mm | 35 mm |
| 130-225-90 | Medium  | 95 mm  | 35 mm |



Weisman Speculum 90° (Left Open)  
With Smoke Tube

| Cat #      | Size    | Length | Width |
|------------|---------|--------|-------|
| 130-226-90 | X-Large | 135 mm | 40 mm |
| 130-220-90 | Large   | 115 mm | 35 mm |
| 130-224-90 | Medium  | 95 mm  | 35 mm |



Cusco Rocket Speculum  
With Smoke Tube

| Cat #   | Size                   | Length | Width |
|---------|------------------------|--------|-------|
| 130-276 | Virgin                 | 86 mm  | 16 mm |
| 130-278 | Small                  | 95 mm  | 27 mm |
| 130-280 | Small Long             | 107 mm | 27 mm |
| 130-282 | Small Extra Long       | 125 mm | 27 mm |
| 130-284 | Medium                 | 95 mm  | 30 mm |
| 130-286 | Medium Long            | 110 mm | 30 mm |
| 130-288 | Medium Extra Long      | 125 mm | 30 mm |
| 130-290 | Large                  | 110 mm | 35 mm |
| 130-292 | Large Long             | 120 mm | 35 mm |
| 130-294 | Extra Large Extra Long | 130 mm | 36 mm |



Winterton Speculum  
With Smoke Tube

| Cat #   | Size  | Length | Width |
|---------|-------|--------|-------|
| 130-296 | Large | 125 mm | 28 mm |

# ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY



Sims Speculum (Large)  
130-240



Sims Speculum (Medium)  
130-242



Sims Speculum (Small)  
130-244



Auvard Weighted Speculum  
130-265



US Army Retractor (Set of 2) 8 1/2" (21.5 cm)  
130-370



Farabeuf Retractor (Set of 2) 6" (15 cm)  
130-380