



**Quality First
Develop Harmoniously**

Rev. No.:2026A

EN ISO 13485:2016

CE0197

Ovlašćeni zastupnik i uvoznik:

MEDICALTECH SOLUTIONS DOO
Ćirpanova 15, Novi Sad
Tel. +381 21 2715 767
Mobilni: +381 62 281 640
office@medicaltech.rs
www.medicaltech.rs


CITEC™

LAPAROSCOPIC PRODUCTS

Date of issue: 2026-01-27

Changzhou Intl. Trade & Enterprises Cooperative Co., Ltd(CITEC)

Disposable TPU Retrieval Bag



Features:

- + High strength of TPU material of the bag ensures no leakage.
- + Self-opening pouch design provides a rapid safe and easy removal of specimen.
- + Transparent bag provides excellent visibility of specimen.
- + Convenient delivery introducer contributes to rapid product entry into abdominal cavity.
- + Various of sizes available fit for different surgery requirements.
- + The model with string provides convenience for surgeons.

Code No.	Size of Bag(mm)	Bag Volume (ml)	Units/ Box
CIT-RB-200-S (With String)	120 x 180	200	10
CIT-RB-350-S (With String)	160 x 200	350	10
CIT-RB-480-S (With String)	200 x 210	480	10
CIT-RB-700-S (With String)	210 x 200	700	10
CIT-RB-200	120 x 180	200	10
CIT-RB-350	160 x 200	350	10
CIT-RB-480	200 x 210	480	10
CIT-RB-700	210 x 200	700	10

CE0197

Endobag



With String

Features:

- + High strength TPU material of the bag ensures no leakage.
- + Transparent bag provides a better visibility of the specimen.
- + Preloaded memory wire automatically opens the bag upon deployment.
- + The string helps to close and retrieve the bag easily (EB-200, EB-400).
- + After deployment, the introducer can be withdrawn, free the cannula.
- + New U shape bag design, disperse the air pressure inside the bottom of bag when pulling out the bag from the body.

Code No.	Introducer Diameter(mm)	Size of Bag(mm)	Bag Volume (ml)	Units/Box
BBD-TSG-01	10	95 x 195	200	10
BBD-TSG-02	10	130 x 195	400	10
BBD-TSG-03	10	160 x 200	800	10
EB-200 (with String)	10	95 x 195	200	10
EB-400 (with String)	10	130 x 195	400	10